



# Medical Equipment Repair Associates

Dedicated to service, integrity and professionalism

## SERVICE REPORT

Workorder # 183299

Page 2 of 3

MERA® Office	TSP-VT	CUSTOMER	NAME	UNIVERSITY OF VERMONT	
Address	280 EAST AVENUE		STREET	SUITE/DEPT	
	BURLINGTON, VT 05401		180 COLCHESTER AVE		
Phone	802-656-3255		CITY	STATE	ZIP
Date Taken	07/16/08 <input type="checkbox"/> Complete		BURLINGTON	VT	05405
		PHONE	CONTACT		
		802-656-3519	JANE DUXBURY		

Date(s) on Site 07/18/08

MANUFACTURER & MODEL #	SERIAL NUMBER	ACCOUNT	SHELDON MANUFACTURING, INC.
SHELDON	10058207		800-322-4897
FX28-2		AUTH/P.O.	
Equipment Description:	BILLING <input checked="" type="checkbox"/> Warranty <input type="checkbox"/> Contract <input type="checkbox"/> Non-Warranty <input type="checkbox"/> Billable		

Problem Reported / Found: DISPLAY IS READING LOW TEMP.

Action Taken: REPLACED TEMP SENSOR, POWERED UP UNIT, LET UNIT STABILIZE, DISPLAY STILL READING LOW, SPOKE WITH JUAN AT SHELDON MANUFACTURING JUAN WILL SHIP OUT TEMP BD SECOND DAY.  
 07/23/08 REPLACED TEMP CONTROL BD, DISPLAY DID NOT COME ON, WORKED WITH SHELDON TECH, SHELDON WILL SHIP OUT NEW PART OVER NIGHT.  
 07/24/08 INSTALLED CONTROL BD FUNCTION TESTED PASSED RETURNED TO SERVICE..

PARTS USED				BREAKDOWN					
Part Number	Qty	Description	Price	LABOR	Standard	TRAVEL	Standard	MILEAGE	Standard
					Overtime		Overtime		Zone
					Special		Special		
					Per Diem		Tolls		Parts
					Airfare		Parking		Shipping
Overtime Authorized by:					Car Rental / Gas		Lodging		Tax
					Taxi		Meals		—
<div style="border: 1px solid black; width: 100%; height: 40px; display: flex; justify-content: flex-end; align-items: center; padding-right: 5px;">+</div>									

Customer Signature

Print Customer Name

Technician Name

P.O. BOX 167 • CHETEK, WI 54728 • 800-776-7158 • (715) 925-1298 • FAX (715) 925-1079

MERA® is a registered service mark of Medical Equipment Repair Associates.

QSF-14 MERA Service Report

Revision 05-0616

Retain:\* 3 years